

**This presentation rated R
for REALLY IMPORTANT**

Report Documentation Page				Form Approved OMB No. 0704-0188	
Public reporting burden for the collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to a penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.					
1. REPORT DATE 05 APR 2008		2. REPORT TYPE		3. DATES COVERED 00-00-2008 to 00-00-2008	
4. TITLE AND SUBTITLE Real-time Monitoring of our Warfighters Health State: The Good, The Bad, and The Ugly				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S)				5d. PROJECT NUMBER	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) U.S. Army Research Institute of Environmental Medicine (USARIEM), Natick, MA, 01760-5007				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)	
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION/AVAILABILITY STATEMENT Approved for public release; distribution unlimited					
13. SUPPLEMENTARY NOTES Army Telemedicine Partnerships Series 2008: ?Personal Health Monitoring?, Seattle, WA, 5 Apr 2008					
14. ABSTRACT					
15. SUBJECT TERMS					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT Same as Report (SAR)	18. NUMBER OF PAGES 28	19a. NAME OF RESPONSIBLE PERSON
a. REPORT unclassified	b. ABSTRACT unclassified	c. THIS PAGE unclassified			

Real-time Monitoring of our Warfighters Health State:

Presented by
COL Beau Freund
Commander USARIEM
at ATA 2008
Seattle, WA

R This presentation rated R
for REALLY IMPORTANT




Disclaimer

- The opinions or assertions contained herein are the private views of the author(s) and are not to be construed as official or reflecting the views of the Army or the Department of Defense.




Disclaimer 2

- The idea of monitoring a Warfighter's health state is not new.



Circa 1956



Circa 1959



Background: Why real-time monitoring?

- Situational awareness is limited to the immediate presence of a medic
- No tool to monitor individual/unit health readiness or triage casualties at a distance
- No remote awareness of casualty vital signs
 - Medics attending to Soldiers KIA
 - endangers medics
 - limits resources



Background: Why real-time monitoring?



- Thermal casualties
 - reduce unit readiness
 - can be fatal
- Dehydration
 - ↑ risk thermal casualty
 - ↓ physical & cognitive performance
- Limited Sleep
 - ↓ cognitive performance
 - ↑ preventable accidents

Paramount Considerations

- Minimal Power, Weight, and Cube
- Open Architecture to Facilitate Upgrades
- Block/Spiral Development
- Keep it Small and Simple
- Minimize Logistics Trail
- Economically Reproducible



Decision Assist Tool Complexity

- **No medic eyes on**
 - **Hard to test Vital Sign Algorithms**
 - **False positives must be very low**
 - i.e. Don't say someone is dead when they are alive and can be helped
 - **Human device interface is complex**
 - **Artifact**
-
- | | | |
|---------------------------|------------------------|---------------------------------|
| ■ Hardware | ■ Activity | ■ Goop & Goo |
| ■ Is it connected? | ■ Movement / placement | ■ Grease / Camouflage |
| ■ How is the battery? | ■ Talking | ■ Repellents (DEET, promethryn) |
| ■ Did it slip? | ■ Vehicle movement | ■ Sugars |
| ■ Proper tension? | ■ Gun fire | ■ Washing / Cleaning |
| ■ Was it shot off? | ■ Run with load | ■ Chlorine |
| ■ Was it properly put on? | ■ Pack shifting | ■ Detergent |



Concept

Medic



Physiological Monitor

Health Status





The GOOD

Technologies & a solution framework have been greatly advanced. Health state monitoring is no longer science fiction.

Maturing Sensors
WPSM Design Principles
WPSM Design Framework
Certification and Lab Testing
Medical Monitoring



Maturing Sensors

- **Sensor Development**

- **Foster Miller**

- *Added ECG and Respiration Waveform Capabilities*

- **Hidalgo System**

- *Added Core Temperature Capabilities*

- **Quasar**

- *Fielding a initial system to Aberdeen Test Center (Jan/Feb 2008)*

- **Zephyr**

- *Commercially available system with ECG and Resp. Waveform*



Foster Miller



Hidalgo Equivital



Quasar



Vivometrics
Lifeshirt 300S



Zephyr

**The
GOOD**


Real-Time Monitoring Design Principles

- Modular
- Open
- Extensible
- Personalizable
- Spiral Development
- Initial Capability (IC) “Locked”



**The
GOOD**

Real-Time Monitoring Design Framework

Current	Sensors	Network	Info. Management	Customers
<i>Real-Time Monitoring</i>	VSDS BIDS FIM Core Pill Temp Patch	Hardwire Mini Mitter Vital Sense Network Blue Tooth (FFW)	HEALTH HUB State Classifiers: Vital signs, Thermal, Hydration Sleep Performance Models	Land Warrior Future Force Warrior
<i>P3I</i>  <i>Future</i>	<ul style="list-style-type: none"> • Non-Contact ECG • Cardiac Output • Non Invasive Core Temp. • Cardiac Output 	<ul style="list-style-type: none"> • Inductive • Zigbee • Stealthy • Through Body • Ultra-wide Band 	<ul style="list-style-type: none"> • Movement Classifiers • Metabolic Rate • Individual Models • BIDS Wound Location • BIDS Wound Severity • Triage 	<ul style="list-style-type: none"> • Early Insertion Efforts • SOF • Rangers • CST-WMD • Training Units



The
GOOD

Real-Time Monitoring: Initial Capability System



BMIS-T

Medic Display



Hub

911 Button

Runs Health state algorithms

Interfaces with the PIC/EIC

FFW Link

→ GPS, Environmental

← Health State Metrics



Vital Sign Detection System (VSDS)

ECG, Respiration, Ambulation, Body Orientation,
Skin Temp.



Core Temp

Thermal State



Fluid Intake Monitor

Hydration State



Sleep Watch

Cognitive/sleep Performance State



The
GOOD

Certification and Lab Testing



■ FDA 510K Certified (November 2006):

- Heart Rate Monitor
- Respiration Effort Monitor
- Welfare Index

■ Lab Validation Testing

- February 2006
 - HR +/- 2 bpm
 - RR +/- 3 bpm (for low activity)



■ Human Factors Testing

- 2004 early form factor testing
- Feb & June (2006) Feb (2007)
 - Devices worn during field exercises for 8 hours – 72 Hours
 - Iterative design process to fix comfort and durability issues



**The
GOOD**

1st CST-WMD Real-Time Medical Monitoring

Methods

1st CST-WMD Enclosed Space Training Event (June 7 & 8, 2007)

- 2 Day Event
- Monitor Environment
- Monitor Physiology
- Display Physiology in Real Time for Use by Physician Assistant





The
GOOD







Medical Monitoring Telemetry System – In Action

Results – Physiology, Real Time Display

CST Medical Telemetry Monitor

Team Detail Diagnostics

 **1st CST - WMD** 
Medical Telemetry Monitor

State	ID	Name	HR		Br/min		Temp.		Orientation	Last Update
✓	6	3060005	131	→	27	→	35.3	→		11:24:32
✓	2	3060001	127	↑	39	→	30.4	→		11:24:22
✓	3	1150001	121	→	32	↑	36.1	→		11:24:29
✓	4	3060004	90	→	36	↑	31.59	→		11:24:32
✓	5	3060008	96	→	28	↓	30.09	→		11:24:19
✓	1	3060012	134	→	25	↓	34.7	→		11:24:25

OK | Current Time: 11:24:35

C R S





The BAD

Unique challenges make the creation and fielding of a product difficult.

Changing Objectives
Changing Requirements
Changing Configurations



Changing Objectives

- **2004 Comprehensive Physiological Status Monitoring System – Integrated into Objective Force Warrior (OFW)**
 - Thermal State, Hydration State, Cognitive State, Vital Sign State, Ballistic Impact Detection
 - Technology Down Select
 - Human Factors form fit design to meet OFW ensemble
 - Integrated algorithms and network
- **2005 Stand Alone WPSM system for Future Force Warrior (FFW)**
 - Intellectual Property and Funding Issue to Integrate Network
 - FDA Issue to integrate algorithms
 - Hydration State requirement dropped
- **2006 Stand Alone Vital Sign Detection System (VSDS)**
 - Change from WPSM network to standard Bluetooth Network
 - Change to new interface specifications
- **2007 FFW Program Ends - Now Integrate to Land Warrior**
 - Hardwired interface required because of security issues on the battlefield
 - Only leaders will have radio systems

Two original ATO's schedules were synchronized.

The FFW ATO was extended but the WPSM ATO was not.

Thus WPSM and FFW schedules became out of sync.



Changing Requirements

VSDS was originally specified with an integrated Ballistic Impact Detection Technology

- Final specification was simplified
 - Only Vital Signs Detection was wanted/affordable

The communication link for health status information was planned as the individual radio

- Land Warrior: current decision is to have radio only on the team leaders




**The
BAD**

Changing Configuration Timeline

Initial Design	Future Force Warrior (FFW) Design	Land Warrior (LW) Design
Multiple sensors	One sensor	One sensor
Wireless sensor network	Bluetooth wireless sensor connection	Wired sensor connection
No Soldier computer	Soldier computer	Soldier computer and radio only on leaders
No Soldier radio	Soldier radio	Interceptor body armor with no standoff
No body armor	Integrated load carriage and body armor with standoff	Side sensor location
PDA display	Front chest sensor location	





The UGLY


The DOD acquisition process is complex. Redirecting early component work to enable a good fieldable system requires time, effort, and coordination.

Deliberate Capabilities Development process
FFW System of Systems




The
UGLY

FFW “System of Systems” (SOS)



The Deliberate Capabilities Development Process (JCIDS)

Integrated Defense Acquisition, Technology, & Logistics Life Cycle Management Framework



Summary

- Despite current challenges of fielding integrated solutions the future is very bright.
- Significant advances in technologies and supporting firmware/software will continue to be made
- Academia, industry and the Defense Department will share in the further development and exploration of these technologies.
- The potential of future applications to home health care, safety monitoring, professional & recreational athletes, research and the Department of Defense are tremendously exciting.



Summary

- Lab on a chip
- The individual “dashboard”
- Personalized medicine
- Remote “office calls”
- Affordable, integrated solutions
- Performance optimization based on individual’s physiology
- Previously unattainable real-time data sets.



Summary

50 years ago Starship Troopers was written as science fiction

What of today's science fiction will be possible in 2059?



CAST

Dr. Reed Hoyt
MAJ Bill Latzka
Mr. Mark Buller
Dr. Jacques Reifman
COL Karl Friedl
Dr. Jack Obusek
COL Bob Vandre
Mr. Tony Karis
Mr. Steve Mullen
Mr. Bill Tharion
Mr. John Ames
Dr. Vic Convertino
Dr. Miyo Yokota
Mr. Steve Van Albert
Dr. Fred Pierce
MAJ Kaushek

